

Johnson[®]

TEST PAPERS



WWW.ABMOLI.COM
09-7434264 | 072-2127007



pH Strips (non-bleeding)

The indicator pads on these environmentally friendly strips are prepared as a non bleed system therefore the resultant colour change remains far longer and readable until the pad is dry. Strips with the non-bleed system provide precise pH values as the different colours do not mix at the point of testing. Each strip is long enough to protect the user from the test solution as the test pads are at the extreme end of the strip. For accurate pH readings, these strips use up to 4 different indicator pads and the colours on the enclosed colour chart match the colour and position of each pad on the strip. This allows these strips to provide a rapid method of measuring the pH of a solution while producing high quality results each time.

Universal Indicator Strips are also available as single pad strips where it is not essential to have the high level of accuracy achieved using the 4 or 3 indicator pad strips. These strips still provide the benefits of using pH indicator strips along with immediate colour changes and a non bleed system to identify the resultant pH.



| Product | Description | Presentation | Product Code |
|--|---|---------------|--------------|
| Universal Indicator Strips pH 0 - 14, 4 pad | pH 0 - 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9 - 10 - 11 - 12 - 13 - 14 | 100 strips/pk | 140.4 |
| Universal Indicator Strips pH 0 -14, 3 pad | pH 0 - 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9 - 10 - 11 - 12 - 13 - 14 | 100 strips/pk | 101.3C |
| Universal Indicator Strips pH 0-14, 1 pad | pH 0 - 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9 - 10 - 12 - 14 | 100 strips/pk | 046.3 |
| Universal Indicator Strips pH 0-11, 1 pad | pH 0 - 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9 - 10 - 11 | 100 strips/pk | 004.3 |
| pH Indicator Strips pH 11.5 - 14.0 | pH 11.5 - 12.0 - 12.5 - 13.0 - 13.5 - 14.0 | 100 strips/pk | 111.3C |
| pH Indicator Strips pH 10.0 - 13.0 | pH 10.0 - 10.5 - 11.0 - 11.5 - 12.0 - 12.5 - 13.0 | 100 strips/pk | 110.3C |
| pH Indicator Strips pH 7.0 - 14.0 | pH 7.0 - 7.5 - 8.0 - 8.5 - 9.0 - 9.5 - 10.0 - 10.5 - 11.0 - 11.5 - 12.0 - 12.5 - 13.0 - 13.5 - 14.0 | 100 strips/pk | 109.3C |
| pH Indicator Strips pH 7.0 - 10.0 | pH 7.0 - 7.5 - 8.0 - 8.5 - 9.0 - 9.5 - 10.0 | 100 strips/pk | 127.2C |
| pH Indicator Strips pH 5.0 - 9.0 | pH 5.0 - 5.5 - 6.0 - 6.5 - 7.0 - 7.5 - 8.0 - 8.5 - 9.0 | 100 strips/pk | 126.2C |
| pH Indicator Strips pH 4.0 - 10.0 | pH 4.0 - 4.5 - 5.0 - 5.5 - 6.0 - 6.5 - 7.0 - 7.5 - 8.0 - 8.5 - 9.0 - 9.5 - 10.0 | 100 strips/pk | 106.3C |
| pH Indicator Strips pH 4.0 - 7.5 | pH 4.0 - 4.3 - 4.6 - 4.9 - 5.2 - 5.5 - 5.8 - 6.1 - 6.4 - 6.7 - 7.0 - 7.5 | 100 strips/pk | 125.2C |
| pH Indicator Strips pH 3.0 - 6.0 | pH 3.0 - 3.5 - 4.0 - 4.5 - 5.0 - 5.5 - 6.0 | 100 strips/pk | 144.1C |
| pH Indicator Strips pH 2.5 - 4.5 | pH 2.5 - 3.0 - 3.3 - 3.6 - 3.9 - 4.2 - 4.5 | 100 strips/pk | 114.2C |
| pH Indicator Strips pH 2.0 - 9.0 | pH 2.0 - 2.5 - 3.0 - 3.5 - 4.0 - 4.5 - 5.0 - 5.5 - 6.0 - 6.5 - 7.0 - 7.5 - 8.0 - 8.5 - 9.0 | 100 strips/pk | 113.3C |
| pH Indicator Strips pH 1.0 - 5.0 | pH 1.0 - 1.5 - 2.0 - 2.5 - 3.0 - 3.5 - 4.0 - 4.5 - 5.0 | 100 strips/pk | 148.1C |
| pH Indicator Strips pH 0 - 6.0 | pH 0 - 0.5 - 1.0 - 1.5 - 2.0 - 2.5 - 3.0 - 3.5 - 4.0 - 4.5 - 5.0 - 5.5 - 6.0 | 100 strips/pk | 103.3C |
| pH Indicator Strips pH 0 - 3.0 | pH 0 - 0.5 - 1.0 - 2.0 - 3.0 | 100 strips/pk | 143.1C |
| pH Indicator Strips pH 0 - 2.5 | pH 0 - 0.5 - 1.0 - 1.5 - 2.0 - 2.5 | 100 strips/pk | 122.2C |
| pH Indicator Strips pH 0 - 1.5 | pH 0 - 0.3 - 0.5 - 0.8 - 1.0 - 1.5 | 100 strips/pk | 121.2C |

